


CHARACTERISTICS

- Protects the patient during oxygen therapy
- The Hydrocolloid material is gentle to the skin
- Contact between the nasal cannula and the skin is minimised
- Holes are significantly larger than the cannula thereby providing a clear flow path

- A hook and loop system securely holds the nasal cannula while enabling its position to be easily adjusted.
- Observation of the condition of the nasal nares and repositioning is easy.
- Reduces irritation caused by shear motion.



Conlett

U25/17 Cairns St
Loganholme
QLD 4129
Australia

Tel: (+61) 7 3209 9957
Email: contact@conlett.com

Description	ORDER#	QTY/BOX
High Flow Neo-Guard Small	HFNG-S	50
High Flow Neo-Guard Medium	HFNG-M	50
High Flow Neo-Guard Large	HFNG-L	50